KeyDenture® Base

(FR)

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KeyDenture® Base

(EN)

Product Description: KeyPrint® KeyDenture® Base is designed for additive manufacturing in vat polymerization 3D printers utilizing wavelengths between 385nm-405nm. KeyDenture® Base is intended to be used within a computer-aided design and manufacturing (CAD/CAM) digital dentistry system that includes a 3D scanner, design software, 3D printer, and post-cure unit for the fabrication of oral denture appliances for fully and partially edentulous patients. For any components that are used in conjunction with KeyDenture® Base, the user should review all applicable product labeling including instructions for Use, user manuals, and other associated labeling. Strict adherence to all labeling is critical in assuring a safe and effective printed appliance.

Indications for Use: KeyPrint KeyDenture® Base is a photopolymer resin indicated for the 3D fabrication of denture bases, including removable full and partial dentures, in dental laboratories. Fabrication requires a computer-aided design and manufacturing (CAD/CAM) system and 3D printer.

Contraindications: Contains acrylate monomers and oligomers which, although rare, may cause an allergic reaction in individuals sensitive to acrylic containing products.

Adverse Reactions: Allergic skin reactions may occur in susceptible individuals. The resin may cause respiratory irritation.

CAUTION: US Federal law restricts this device to sale by, or on the order of a dental professional.

- Warnings & Precautions:**
- Follow all compatible printer settings and curing times for safe and effective results. Refer to the QR codes below.
 - Do not use any devices or components that are not KeyPrint compatible or deemed acceptable by Keystone Industries. Refer to the QR codes below.
 - Review the product Safety Data Sheet (SDS) prior to use.
 - As per the SDS, wear proper personal protective equipment when handling KeyPrint resins and uncured printed parts.
 - Agitate the bottle 1 hour prior to use for a minimum of 5 minutes to ensure all pigments are dispersed within the resin and to avoid bubbles during printing.
 - When pouring the resin, be careful not to splash.
 - Store in a cool, dry place 15°C-30°C (59°F-86°F) and away from light. Ensure that the bottle is capped while not in use.
 - In the unlikely event of a print failure during printing, filter the liquid resin through a mesh screen with pore sizing <200 microns. It is good practice to filter the resin in a vat periodically to prevent print failures.
 - Clean out the printer platform and vat prior to using a different batch of resin. Keystone recommends designating a vat that is specific for each KeyDenture® shade. DO NOT mix different batches of the same product.
 - Do not pour any remaining print resin material back into the bottle after printing.

Compatible Equipment: To ensure the biocompatibility of the final device, Keystone Industries collaborates with printer manufacturers to provide compatible printer and post-cure settings. Refer to the QR Codes for a full list of compatible workflows.

- Directions for Use:**
- Ensure the resin is tempered to ambient temperature (20-25°C/68-77°F) prior to printing.
 - Agitate the bottle 1 hour prior to use for a minimum of 5 minutes to ensure all pigments are dispersed within the resin and to avoid bubbles during printing.
 - Only use KeyDenture® compatible printers, settings, and post cure units to ensure the printed product is safe and accurate for use. Refer to the QR codes above.
 - For printer settings see instructions for use (IFU) of 3D printer. Follow the manufacturer's IFU of the printer.
 - Once your design is completed per CAD software manufacturer's directions for use, import the design into the CAM software unique to the printer manufacturer.
 - Clean the printed part per the instructions below.

Cleaning Instructions for Printed Parts: Printed parts should be cleaned with Isopropanol alcohol (IPA) (>97% purity) within approximately eight (8) hours from completion of the print using the following steps:

* Keystone discourages the use of denatured alcohol or ethanol for cleaning as they may diminish or degrade the quality of the finished parts.

- Using a soft brush or agitation, wash the printed parts in clean IPA (>97% purity) bath for 2 minutes.
- Remove printed parts from initial bath, blow dry individual components using compressed air, place in another clean bath of IPA (>97% purity), and wash for an additional 2 minutes.
- Remove printed parts from the bath and blow dry using compressed air.
- Ensure all residual resin is removed from the parts, as residual resin will become shiny and highly noticeable. If resin is still present, use a soft brush to clean that isolated area.
- Allow the IPA to evaporate before post curing the printed parts. A minimum wait time of 30 minutes should be used.

Note: Printed parts should not exceed 5 minutes of IPA exposure.

- Directions for post cure of printed part(s):**
- Refer to the QR code above for all compatible workflows.
 - For optimal device quality, allow parts to cool completely before handling.

At this stage the medical device is cured and safe with respect to residual monomers.

Patient Cleaning Instructions: This medical device is a single-patient, customized, multi-use oral appliance that should be cleaned between usages. The patient should clean the appliance with soap and warm water, or any over-the-counter mild cleaning agents indicated for oral devices.

Disposal Considerations: Unused resin materials may be considered hazardous waste, dispose of in accordance with all local, regional, national, and international regulations.

KeyDenture® Base

Compatible Equipment
Compatibele apparatuur
Équipement compatible
Kompatible Geräte
Apparecchiature compatibili
Equipamento compatível
Equipo compatible

In case of Emergency:
INFOTRAC: 800-535-5053
Customer# 76987

■ Keystone Industries
52 West King St.
Myerstown, PA 17067 USA
800.333.3131
keystoneindustries.com

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